Ruhvenile® Biomedical
Product datasheet

ViroStabl™ UTM

Universal Transport Medium for Viral Culture and Molecular Testing

Keep pathogens viable for culture and intact for molecular analysis, with 5-day room-temperature transport.

Indigenous patented technology | Discovered and developed in India | DSIR recognised | ISO 9001:2015 certified | CE (IVD) compliant



120 hours	1 year	12–18 months	2.0 mL
transport at 2–25°C	at -80°C	shelf life at room temperature	medium in a 5.0 mL tube

Key features

- Up to 120 hours of transport at 2–25°C without a cold chain.
- Supports viral culture, genetic testing and molecular analysis from the same specimen.
- Covers SARS-CoV-2, influenza A/B, RSV, rhinovirus, chlamydiae, mycoplasma and ureaplasma.

Specifications

Parameter	Value
Presentation	2.0 mL medium in a 5.0 mL sterile collection tube; kit with nasal swab available
Specimen stability	Room temperature up to 5 days; -80°C 1 year
Medium shelf life	12–18 months at room temperature, away from direct sunlight
Sample types	Nasal, throat and other validated swab specimens

Catalogue codes and pack sizes

Code	Pack
RVN-VS-2.0-00	2.0 mL in 5.0 mL tube with swabs

Best for

- Labs that need viral culture and molecular testing from the same swab
- Chlamydia, mycoplasma and ureaplasma work
- Respiratory surveillance: SARS-CoV-2, influenza A/B, RSV and rhinovirus

Frequently asked questions

How long can swabs travel in ViroStabl™?

Up to 120 hours (5 days) at 2–25°C.

Can the same swab be cultured and tested by PCR?

Yes. It supports viral culture and molecular analysis from the same specimen.

Which organisms does ViroStabl™ support?

SARS-CoV-2, influenza A/B, RSV, rhinovirus, chlamydiae, mycoplasma and ureaplasma.

What is the shelf life of ViroStabl™?

12–18 months at room temperature, away from direct sunlight.

Is a swab kit available?

Yes, a kit with a nasal swab is available.

Order or customise

Request a quote, a custom pack size or bulk quantities. Dispatch within 3–5 working days.

Ruhvenile® Biomedical OPC Pvt. Ltd.
Mobile and WhatsApp: +91 858 703 3367 | Landline: +91 11 3501 7209
info@ruhvenile.com | www.ruhvenile.com/virostabl.html

Specifications come from Ruhvenile® product information and may be updated. Confirm current specifications and regulatory status with Ruhvenile® before ordering.